



Excellence in Pathology and Laboratory Services

Test Update

From Your Laboratory Service Provider, PathGroup Labs

Methodology change for Vitamin D 25-Hydroxy Testing

Samuel Smith, MD Lab Director / Jeff Johnson, Clinical Operations Manager

Overview and Clinical Utility:

Effective September 1, 2010 PathGroup Labs will begin performing Vitamin D 25 OH using chemiluminescence technology. Test code ordering will remain the same. Report format will also remain the same. Turnaround time will decrease with the new methodology, allowing for faster turnaround time for results.

Test Methodology:

This chemiluminescence technology on the IDS iSYS Multi-Discipline Automated Analyzer is intended for the quantitative determination of 25-hydroxyvitamin D (25OHD) in human serum. In conjunction with clinical findings, this test is indicated to assist the clinician in the assessment of vitamin D sufficiency in an adult population.

Ordering:

Test Code	Test Name:	CPT:	LOINC:
VD25H	Vitamin D 25-Hydroxy	82306	1989-3

Reference interval for 25 OHD:

New Reference Range:

Methodology: Chemiluminescence	Methodology: EIA
30 - 80	30 - 80
Units: ng/mL	Units: ng/mL

Interpretation of Vitamin D 25 OH levels:

Level	Range (ng/mL)
Deficient	< 10
Insufficient	10 - 29
Sufficient (Optimal)	30 - 100
Potential Intoxication	> 100

Specimen Collection and Storage:

- 1) Collect specimen in a serum separator tube.
 - a) Allow serum tube to clot for 30 minutes.
 - b) Spin serum separator tube 10-15 minutes @ 3500 rpm.
- 2) Serum specimens are stable refrigerated for 3 days or frozen for 30 days.

Test Performed:

Monday - Friday

Turnaround time:

24 hours

For further questions, please contact Client Services at 615-562-9300 or 1-888-474-5227.