

## TEST UPDATE SERUM IODINE

---

Mary M. Mayo, PhD, Medical Director

### OVERVIEW AND CLINICAL UTILITY:

**Effective April 27, 2015**, PathGroup will implement updated specimen collection requirements for Serum Iodine testing. As a result, plasma is no longer an acceptable specimen of choice for Serum Iodine testing.

### ORDERING:

TEST CODE	TEST NAME	CPT CODE
IODSE	Iodine, Serum	83789

### SPECIMEN COLLECTION AND STORAGE:

#### IF ORDERING:

- Royal Blue, no additive
- Separate serum from cells within 2 hours of collection
- Aliquot 2ml of serum
- Transport Refrigerated

**For further questions, please contact Client Services at 615-562-9300 or 1-888-474-5227.**