

Test Update: Specimen Requirement Changes

November 6, 2019

Effective November 17, 2019, the test codes below will have changes to specimen containers and/or transportation temperatures. If specimen collection containers for a test have changed, you may need to order new supplies from PathGroup prior to the effective date. For additional information, please refer to the PathGroup Directory of Services at <http://www.pathgroup.com/clinical/testing-menu/>.

Test Code	Test Name	Specimen Container		Transportation Temperature	
		Old	New	Old	New
17HPS	17-Hydroxypregnenolone Quant, MS/MS, Serum	Serum Separator Tube, Pour Off	Red-Top Tube, Pour Off	Frozen	No Change
FETYL	Fentanyl and Metabolite, Serum or Plasma, Quantitative	Red-Top Tube, Pour Off	EDTA (Lavender), Whole Blood	Refrigerated	No Change
OHPC2	17-Hydroxypregnenolone Quant, MS/MS, Serum	Serum Separator Tube, Pour Off	Red-Top Tube, Pour Off	Frozen	No Change
HSPE	Hypersensitivity Pneumonitis (Farmer's Lung) Extended Panel	Serum Separator Tube, Pour Off	No Change	Refrigerated	Frozen

For further questions, please contact Client Services at 615-562-9300 or 1-888-474-5227.