

Test Update: Specimen Requirement Changes

November 26, 2019

Effective December 8, 2019, the test codes below will have changes to specimen containers and/or transportation temperatures. If specimen collection containers for a test have changed, you may need to order new supplies from PathGroup prior to the effective date. For additional information, please refer to the PathGroup Directory of Services at <http://www.pathgroup.com/clinical/testing-menu/>.

Test Code	Test Name	Specimen Container		Transportation Temperature	
		Old	New	Old	New
C1EIF	C-1-Esterase Inhib. Functional	Serum Separator Tube	Red-Top Tube, Pour Off	Frozen	No Change
C1EIP	C1 Esterase Inhibitor Panel, Antigen and Functional	Serum Separator Tube, Pour Off	Red-Top Tube, Pour Off	Frozen	No Change
HEMUR	Hemosiderin, Urine	Sterile Container	Sterile Container	Frozen	Refrigerated
INSFB	Insulin, Free and Total	Serum Separator Tube	Serum Separator Tubes x2	Frozen	Refrigerated
KSUSP	Urine Supersaturation Profile	Sterile Container	24 Hour Aliquot x2 (min 30 mL each)	Frozen	Ambient

For further questions, please contact Client Services at 615-562-9300 or 1-888-474-5227.