

## Test Update: Specimen Requirement Changes

December 23, 2019

The test code has had changes to specimen containers and/or transportation temperatures. If specimen collection containers for a test have changed, you may need to order new supplies from PathGroup. For additional information, please refer to the PathGroup Directory of Services at <http://www.pathgroup.com/clinical/testing-menu/>.

Test Code	Test Name	Specimen Container		Transportation Temperature	
		Old	New	Old	New
PH	pH, Urine, 24-hr Urine Collection	Urine Chemistry Red-Top Tube	UA Preservative Tube	Refrigerated	No Change

**For further questions, please contact Client Services at 615-562-9300 or 1-888-474-5227.**