
Test Update: Specimen Requirement Change

May 5, 2020

Effective immediately, the **Tick Identification** test (**Test Code: TICKID**) should be collected and transported in saline or ethanol in a clean container. Transportation liquids like formalin will no longer be accepted for testing.

| Test Code | Test Name | Specimen Container | |
|-----------|--|--------------------|-------------------|
| | | Old | New |
| TICKID | Tick (and other Arthropods) Identification | Formalin | Saline or Ethanol |

**For further questions, please contact Client Services
at 1-888-4PG-LABS (1-888-474-5227).**